



- 8) McCoy MEGA Blades - FIBEROPTIC**
- 9) McCoy Macintosh Blades - Standard**
- 10) McCoy Macintosh Blades - FIBEROPTIC**
- 11) Miller Optima CLX Blades - FIBEROPTIC**
- 12) Other Blades**
- 13) Disposable Blades - Standard**
- 14) Disposable Blades - FIBEROPTIC**

FlexTip 70°



McCOY - MEGA Blade

Fiber Optic illumination



Specification:

The McCoy Is A Macintosh Laryngoscope Blade With A Tip Which Is Adjustable Through 70°. A Lever Controls The Tip Angle During Incubation To Lift The Epiglottis, Giving A Clear View Of The Vocal Cords. Contact With The Upper Row Of Teeth And Unnecessary Pressure On The Tongue And Soft Tissue Are Avoided.

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Inserted Fiber Optic Technology For Advanced Performance Up To 6500 Individual Micro-fibers For Improved Light Transmission And Longer Life.
- ▶ These Blades Need No Maintenance — Ever.
- ▶ Built With An Integral, Extra-large Fiber Optic Bundle Providing Cooler, White Light.
- ▶ Smooth Shape Avoids Sharp Edges And Cavities For Easy Cleaning And Sterilizing..
- ▶ Integral Design Enables Superior Cleaning By Elimination Of Cracks And Crevices.
- ▶ Both Ends Of The Blades Are Protected For Life.
- ▶ Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-106-04	Child	2	102 MM
75-106-05	Medium Adult	3	130 MM
75-106-06	Large Adult	4	155 MM

Fo MCCOY BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 2 Blades Size 2 and Adult Handle	75-420-53
Set of 2 Blades Size 3 and Adult Handle	75-420-54
Set of 2 Blades Size2, 3 and Adult Handle	75-420-55
Set of 4 Blades Size2, 3 4 and Child Handle	75-420-56



FlexTip 70°



**McCoy - Blade
Standard illumination**

Specification:

The McCoy is a Macintosh Laryngoscope Blade with a tip which is adjustable through 70°. A lever controls the tip angle during intubation to lift the epiglottis, giving a clear view of the vocal cords. Contact with the upper row of teeth and unnecessary pressure on the tongue and soft tissue are avoided.

- ▶ Manufactured from Aisi 304 Stainless Steel.
- ▶ All blades compatible with ISO 7376 standard.
- ▶ Available with LED and standard bulb illumination.
- ▶ Innovative design for simple intubation and good visibility.
- ▶ Autoclave able.
- ▶ Bulbs 2.5v Halogen also 3.5v LED with 20000 Lux.
- ▶ These blades need no maintenance — ever but if problem with bulb then it can be replaced.
- ▶ Quick-to-remove lamp at the front of the blade.
- ▶ The rounded edges of the blade reduce the risk of injury.
- ▶ Require conventional handle.

Size Measurement

Code	Blades	B/No	Overall length
75-107-04	Child	2	102 MM
75-107-05	Medium Adult	3	130 MM
75-107-06	Large Adult	4	155 MM

FO MCCOY BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 2 Blades Size 2 and Adult Handle	75-430-53
Set of 2 Blades Size 3 and Adult Handle	75-430-54
Set of 2 Blades Size 2, 3 and Adult Handle	75-430-55
Set of 4 Blades Size 2, 3 4 and Child Handle	75-430-56



FlexTip 70°



**Also
FIX Tube**

McCoy - Blade

Fiber Optic illumination



Specification:

The McCoy Is A Macintosh Laryngoscope Blade With A Tip Which Is Adjustable Through 70°. A Lever Controls The Tip Angle During Incubation To Lift The Epiglottis, Giving A Clear View Of The Vocal Cords. Contact With The Upper Row Of Teeth And Unnecessary Pressure On The Tongue And Soft Tissue Are Avoided.

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Available In 3mm, 4mm, 5mm Fiber Optic Illumination
- ▶ Inserted Fiber Optic Technology For Advanced Performance
- ▶ Individual Micro-fibers For Improved Light Transmission.
- ▶ Smooth Shape Avoids Sharp Edges And Cavities For Easy Cleaning And Sterilizing.
- ▶ Extra-strong Fiber Bundles To Ensure The Best Possible Light Intensity And Durability.
- ▶ The Clear Design Of The Macintosh F.O. Blades Simplifies Cleaning And Autoclaving
- ▶ Fix Tube - Fiber Optic Pipe Is Non-removable It Cant Be Remove.
- ▶ Also Available In Standard Fiber Optic Illumination (removable pipe). Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-105-04	Child	2	105 MM
75-105-05	Medium Adult	3	130 MM
75-105-06	Large Adult	4	155 MM

Fo MCCOY BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 2 Blades Size 2 and Adult Handle	75-440-53
Set of 2 Blades Size 3 and Adult Handle	75-440-54
Set of 2 Blades Size2, 3 and Adult Handle	75-440-55
Set of 4 Blades Size2, 3 4 and Child Handle	75-440-56



Disposable Blades Standard illumination

Specification:

- ▶ Single-use, sterile disposable blades help eliminate the time, reprocessing steps and costs related to disinfecting reusable blades.
- ▶ Packaged individually in a clear pouch making it easy to identify blade size and type.
- ▶ Reduced risk of cross-contamination
- ▶ Stainless steel body
- ▶ Macintosh sizes 0-4
- ▶ Miller sizes 00-4
- ▶ Interchangeable with all standard handles



Miller Size Measurement

Code	Blades	B/No	Overall length
75-200-93	Premature	00	68 MM
75-200-94	Neonate	0	80 MM
75-200-95	Infant	1	95 MM
75-200-96	Child	2	105 MM
75-200-97	Medium Adult	3	135 MM
75-200-98	Large Adult	4	155 MM

Mac Size Measurement

Code	Blades	B/No	Overall length
75-200-99	Neonate	0	80 MM
75-200-00	Infant	1	95 MM
75-200-01	Child	2	105 MM
75-200-02	Medium Adult	3	135 MM
75-200-03	Large Adult	4	155 MM

MACINTOSH CONVENTIONAL SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	75-450-57
Set of 3 Blades Size 1, 2, 3 and Adult Handle	75-450-58
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	75-450-59
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	75-450-60

Disposable Blades Standard illumination

Specification:

- ▶ Single-use, sterile disposable blades help eliminate the time, reprocessing steps and costs related to disinfecting reusable blades.
- ▶ Packaged individually in a clear pouch making it easy to identify blade size and type.
- ▶ Reduced risk of cross-contamination
- ▶ Stainless steel body
- ▶ Macintosh sizes 0-4
- ▶ Miller sizes 00-4
- ▶ Interchangeable with all green-band fiber optic handles



Miller Size Measurement

Code	Blades	B/No	Overall length
75-800-05	Premature	00	68 MM
75-800-06	Neonate	0	80 MM
75-800-07	Infant	1	95 MM
75-800-08	Child	2	105 MM
75-800-09	Medium Adult	3	135 MM
75-800-10	Large Adult	4	155 MM


Mac Size Measurement

Code	Blades	B/No	Overall length
75-800-11	Neonate	0	80 MM
75-800-12	Infant	1	95 MM
75-800-13	Child	2	105 MM
75-800-14	Medium Adult	3	135 MM
75-800-15	Large Adult	4	155 MM

MACINTOSH CONVENTIONAL SETS

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	75-460-57
Set of 3 Blades Size 1, 2, 3 and Adult Handle	75-460-58
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	75-460-59
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	75-460-60



- 15) Standard and F.O. Handles**
- 16) Standard and F.O. Handles - New Pattern**
- 17) Standard and F.O. Handles Rechargeable**
- 18) Charger and Batteries**