

Replacements Bulbs

Bulb	Battery	Code
Xenon	2.5V	75-15
LED	3.0V	75-07
Halogen	2.5V	75-16
Halogen	3.5V	75-160



**Handles
Fiber Optic Illumination**

Specification:

- ▶ Chrome Plated Brass Fiber Optic Handle compatible with all "Green System" fiber optic blades.
- ▶ Require 2 AA size Batteries or 3.5 V Rechargeable Battery.
- ▶ Crafted From High quality Chrome plated Brass.
- ▶ Autoclaveable

Battery Handle Fiberoptic Xenon, Brass Chrome Plated

Handle	Overall length	Battery	Code
Infant	122 MM	1 x AA Size	75-300-17
Stubby	125 MM	1 x C Size	75-300-18
Pediatric	160 MM	2 x AA Size	75-300-19
Medium	160 MM	2x C size	75-300-20
Large	175 MM	2x D size	75-300-21

Battery Handle Fiberoptic LED, Brass Chrome Plated

Handle	Overall length	Battery	Code
Infant	122 MM	1 x AA Size	75-300-22
Stubby	125 MM	1 x C Size	75-300-23
Pediatric	160 MM	2 x AA size	75-300-24
Medium	150 MM	2x C size	75-300-25
Large	175 MM	2x D size	75-300-26

**Handles
Conventional Illumination**

Battery Handle Conventional, Brass Chrome Plated


Handle	Overall length	Battery	Code
Infant	110 MM	1 x AA Size	75-300-27
Stubby	115 MM	1 x C Size	75-300-28
Pediatric	145 MM	2 x AA size	75-300-29
Medium	150 MM	2 x C size	75-300-30
Large	165 MM	2 x D size	75-300-31




Handles
Fiber Optic Illumination

Handle	Overall length	Code
Pediatric	160 MM	75-300-32
Medium	150 MM	75-300-33
Medium With DC	150 MM	75-300-33-1

Rechargeable Handle
Bottom For Desktop Charger



DC Rechargeable Handle
Bottom



Handles
Conventional Illumination

Handle	Overall length	Code
Pediatric	145 MM	75-300-34
Medium	150 MM	75-300-35
Medium With DC	150 MM	75-300-35-1

Handles
Oto & Ophthalmoscope

Handle	Overall length	Code
Aluminum Handle	145 MM	75-300-36
Medium	150 MM	75-300-37

- ▶ Chrome Plated Brass Fiber Optic Handle compatible with all "Green System" fiber optic blades.
- ▶ 3.5 V Rechargeable Battery.
- ▶ Crafted From High quality Chrome plated Brass.
- ▶ Autoclaveable



Desk Charger Twin Port For Mega and Baby Handle

Specification:

- ▶ Manufactured From Plastic Material
- ▶ Big Handles requires 3.5 V size Rechargeable batteries
- ▶ Baby Handles requires 3 V size Rechargeable batteries

Purpose And Features

- ▶ The DUKE is only intended to charge rechargeable battery handles.
- ▶ Power-on indicator. Charging state indicator.
- ▶ When take the handle away the charger indicator
- ▶ Automatically turns off
- ▶ 3. 7V Lithium battery.
- ▶ 2X charging ports.
- ▶ Circuit short, reverse-voltage, over-voltage protection.
- ▶ Low battery indicator.

Use This For Child Handle
Converter for Rechargeable Baby Handle
Charge with Desk Charger



300-38-DUKE
With 9-12 V DC charger
Every Port Provide 4.5 V OUTPUT

Li-ion Batteries



LI-ION BATTERY

For Big Handle 4.0 V
Can Only Charge With Desk
Charger



DC LI-ION BATTERY

For Diagnostics Set & Big Handle
Can Only Charge With Normal
Wall Charger Of 5-6 V
Also Available For Standard And
Fiber Optic Handle



LI-ION BATTERY

For Baby Handle 3.0 V For
Standard And Fiber Optic
Can Only Charge With Desk
Charger

DESCRIPTION	Code
Adult & Otoscope Handle 2.5V	75-300-39
Adult & Otoscope Handle 3.5V	75-300-40

DESCRIPTION	Code
Adult & Otoscope Handle 2.5V	75-300-41
Adult & Otoscope Handle 3.5V	75-300-42

DESCRIPTION	Code
Adult & Small Handle 2.5V	75-300-43
Adult & Small Handle 3.5V	75-300-44



- 19) Oto and Ophthalmoscope Re-usable**
- 20) Otoscope Sets**
- 21) Oto and Ophthalmoscope Sets**
- 22) Dermatoscopes and ENT set**
- 23) Diagnostic Sets**
- 24) Mini Otoscopes**
- 25) Mini Ophthalmoscopes**
- 26) Pocket Style Ophthalmoscopes**
- 27) Pocket Style Dermatoscopes**