

## Otosopes and Ophthalmoscopes

Good Quality and Reusable



75-200-04  
Brass Otoscope Fiber Optic  
Brass Otoscope Conventional



75-200-05  
ENT Set Complete  
Pouch Packing



75-200-06  
Brass Ophthalmoscope

### Specification:

- ▶ Manufactured From Brass Chrome plated
- ▶ 3x Magnifying glass
- ▶ Available in Halogen and LED Illumination
- ▶ Bulbs 2.5V
- ▶ Handles requires 2 x C size batteries
- ▶ Hard Case and Soft Pouch available
- ▶ Fiber Optic Light Source

### Fittings



Pin Fitting



Bayonet Fitting

## Otoscopes And Ophthalmoscopes STANDARD SETS

### BASIC OTOSCOPE SET Bayonet Fitting

Complete With Otoscope, Battery Handle, 1 Set Of 5 Reusable Specula Dia 2.5, 3.5, 4.5, 5.5, 9.0mm In a Pouch

CODE	DESCRIPTION
75-500-72	Completes Set With Spare Lamp
75-500-73	Complete Set Without Spare Lamp
75-500-74	Complete Set With Spare Lamp And Rechargeable Battery
75-500-75	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-76	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger

### BASIC OPHTHALMOSCOPE SET

Complete With OPHTHALMOSCOPE, Battery Handle, In a Pouch

CODE	DESCRIPTION
75-500-77	Completes Set With Spare Lamp
75-500-78	Complete Set Without Spare Lamp
75-500-79	Complete Set With Spare Lamp And Rechargeable Battery
75-500-80	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-81	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger

### STANDARD OTOSCOPE SET Pin Fitting

Complete With Otoscope, Battery Handle, 1 Set Of 5 Reusable Specula Dia 2.5, 3.5, 4.5, 5.5, 9.0mm In a Pouch

CODE	DESCRIPTION
75-500-82	Completes Set With Spare Lamp
75-500-83	Complete Set Without Spare Lamp
75-500-84	Complete Set With Spare Lamp And Rechargeable Battery
75-500-85	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-86	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger

### BASIC NASAL OTOSCOPE SET Pin Fitting

Complete With Otoscope, Battery Handle, 1 Nasal Speculum, In A Pouch  
Also Available In Plastic Box

CODE	DESCRIPTION
75-500-87	Completes Set With Spare Lamp
75-500-88	Complete Set Without Spare Lamp
75-500-89	Complete Set With Spare Lamp And Rechargeable Battery
75-500-90	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-91	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger





## Otosopes And Ophthalmoscopes STANDARD SETS

### BASIC OPHTHALMOSCOPE SET

Complete With Ophthalmoscope, Battery Handle, In A Plastic Box

CODE	DESCRIPTION
75-500-92	Completes Set With Spare Lamp
75-200-06	Complete Set Without Spare Lamp
75-500-94	Complete Set With Spare Lamp And Rechargeable Battery
75-500-95	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-96	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



### STANDARD OTOSCOPE AND OPHTHALMOSCOPE SET

Complete With Otoscope ,ophthalmoscope, Battery Handle, In A Plastic Box

CODE	DESCRIPTION
75-500-97	Completes Set With Spare Lamp
75-500-98	Complete Set Without Spare Lamp
75-500-99	Complete Set With Spare Lamp And Rechargeable Battery
75-500-00	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-01	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



### STANDARD OTOSCOPE SET

Complete With Otoscope, Battery Handle, 1 Set Of 3 Reusable Specula Dia 4.5, 5.5, 9.0mm In Plastic Box

CODE	DESCRIPTION
75-500-02	Completes Set With Spare Lamp
75-200-04	Complete Set Without Spare Lamp
75-500-04	Complete Set With Spare Lamp And Rechargeable Battery
75-500-05	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-06	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



pin Fitting

Bayonet Fitting

## BASIC DERMATOSCOPES SET

Complete With Dermatoscopes, Spare Mirror, In A Plastic Box

CODE	DESCRIPTION
75-500-07	Completes Set With Spare Lamp
75-500-08	Complete Set Without Spare Lamp
75-500-09	Complete Set With Spare Lamp And Rechargeable Battery
75-500-10	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-11	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



## BASIC DERMATOSCOPES SET

Complete With Dermatoscopes, Spare Mirror, In A Pouch

CODE	DESCRIPTION
75-500-12	Completes Set With Spare Lamp
75-500-13	Complete Set Without Spare Lamp
75-500-14	Complete Set With Spare Lamp And Rechargeable Battery
75-500-15	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-16	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



## STANDARD DIAGNOSTIC SET

Complete With Standard Otoscope, Standard Ophthalmoscope, Handle, Nasal speculum, Laryngeal mirrors, Reusable specula, Hard case

CODE	DESCRIPTION
75-500-17	Completes Set With Spare Lamp
75-500-18	Complete Set Without Spare Lamp
75-500-19	Complete Set With Spare Lamp And Rechargeable Battery
75-500-20	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-21	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger





## Diagnostic Set Full Specification

### UNIVERSAL DIAGNOSTIC SET

75-2000-04



- Handle (dry Cell Or Rechargeable) X 1
- Otoscope
- Ophthalmoscope
- Laryngeal Mirrors X2
- Reusable Specula X 5
- Tongue Blade Holder
- Metal Tongue Blade Depressor
- Nasal Speculum
- Ear Tip Connector
- Ear Tips
- Replacement Bulb X 2
- Hard Case X 1
- Also Available In Leather Pouch Just Add -1 After Code

### RIESTER DIAGNOSTIC SET

75-2000-02



- Standard Otoscope X 1
- Standard Ophthalmoscope X 1
- Handle (dry Cell Or Rechargeable) X 1
- Nasal Speculum
- Laryngeal Mirrors 2x
- Reusable Specula X 5
- Replacement Bulb X 2
- Hard Case X 1
- Also Available In Leather Pouch Just Add -1 After Code

### More Details

- Manufactured From High Quality Material
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available

### About Sets

Diagnostic Sets Are The Most Commonly Used Tool Sets During A Clinicians Repertoire To Perform A Comprehensive Physical Test. Diagnostic Sets Facilitate Physicians To Perform Quick Exams And Accurately.



**75-300-04**  
 Mini Otoscope Fiber Optic  
 Complete Set of Disposable Cannulas



**75-300-05**  
 Mini Otoscope Fiber Optic  
 Complete Set of Disposable Cannulas

**Specification:**

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source



A B C D

**Color Options**







75- 500-22



75- 500-20

### Mini Ophthalmoscope

Parker mini Ophthalmoscope with 2.5v LED bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses - 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification fo adopter setting.

Power Supply: 2 x AA Batteries ( not included in set ).  
Packing: Hard plastic case packing.



75-500-23



75-500-24



75- 500-25



### Multi-Action Tongue Depressor

All-Purpose Pocket-Sized Illuminator that is Easy to Use and Maintain.  
Multiple Applications.

Tongue-Depressor Holder can be Used with Disposable Depressors (Included) or with Standard Wooden Depressors. Bright Halogen Bulb with More Intensity Compared with Conventional Bulbs. Refined Plastic Compact Body for Long Lasting Use. Attachment Clip with On/Off Button.

Set Includes: (1) Handle, (1) Head with Integrated Tongue-Depressor Holder (3) Disposable Depressors LED Illuminating Probe. Blue Leather Case. Two "AA" Batteries Required.  
Batteries not Included.

The ophthalmoscope is a tool used in medical to examine the interior of the eye. Studying the inner parts of the eye, such as the optic disk, retinal blood vessels, retina, choroid and macula, can help diagnose diseases of the eye.



## Pocket Ophthalmoscope

75-450

The Pocket size Ophthalmoscope is a lightweight, compact diagnostic instrument that delivers high performance in a slightly smaller package. Quality optics in mini-format is ideal for rounds and house calls. Attachment clip with integrated on/off switch. High-quality plastic handle with shockproof, sturdy, non-slip features. Enhanced LED Technology for perfect illumination and a brilliant image. 4 different apertures. Comprehensive lens range from + 20 D to -20 D.

Power Supply: 2 x AA Batteries ( not included in set ).  
Packing: Hard plastic case packing.



## Pocket Ophthalmoscope

75-500

2.5v LED lamp light for true tissue color and consistent long lasting illumination 16 color coded lenses -25 to +25 for excellent resolution (Red for - Black for +) 4 different apertures (small, large, blue, half moon focus)  
Rubber eye brow rest to prevent scratching of eye-glasses  
illumination lens dial for clear identification of adopter setting.

Power Supply: 2 x AA Batteries ( not included in set ).  
Packing: Hard plastic case packing.





## Dermatoscope

75-550

Modern pocket Dermatoscope in a compact design with high-quality optics. 10X magnification large, high-resolution image. Choice of two skin-friendly and sterilizable contact plates: with scale from 0 - 10mm for exact measurement of pigmented skin lesions, without scale. Simple lamp change on the underside of the instruments head. Enhanced LED Technology for more light. Attachment clip with integrated On/Off switch. High-quality 2 color handle.

Power Supply: 2 x AA Batteries ( not included in set ).  
Packing: Hard plastic case packing.



## Dermatoscope

75-610

Bayonet locking system for easy engagement and disengagement from the handle.

2.5v LED illumination of the examination field. High resolution achromatic lenses system with 10x magnification. Adjustable focusing Eyepiece for a large correction range. Lens without scale, Soft eyepiece protection.

Available in Polish and Sand finish.

Power Supply: 2 x C Batteries ( not included in set ).

Packing: Hard plastic case packing.



- 29) Spare accessories Diagnostic set**
- 30) Bulbs and Spatulas**