

**TONGUE DEPRESSOR**

DESCRIPTION	CODE
Tongue Depressor	75-600-21



**CURVED ILLUMINATOR**

DESCRIPTION	CODE
Bayonet Fitting	75-600-32
Pin Fitting	75-600-33
Re-usable	75-600-34



**METAL TONGUE DEPRESSOR**

DESCRIPTION	CODE
Tongue Depressor	75-600-22



**OPHTHALMOSCOPE HEAD**

DESCRIPTION	CODE
Bayonet Fitting	75-600-35
Pin Fitting	75-600-36
Re-usable	75-600-37



**OTOSCOPE HEAD**

DESCRIPTION	CODE
W/O Ear Spatula	75-600-23
With Ear Spatula	75-600-24
Bayonet Fitting	75-600-25
Screw Fitting	75-600-26
With LED	75-600-27
With Bulb	75-600-28



**NASAL SPECULUM**

DESCRIPTION	CODE
Nasal Speculum	75-600-38



**ADAPTER FOR DISPOSABLE T.D**

DESCRIPTION	CODE
Adapter For Disposable T.D	75-600-29



**SPREADABLE NASAL SPECULUM**

DESCRIPTION	CODE
SPREADABLE NASAL SPECULUM	75-600-39



**METAL ADAPTER FOR SPECULA**

DESCRIPTION	CODE
Metal Adapter For BY Fitting	75-600-30
Metal Adapter For PN Fitting	75-600-31



**LARYNGEAL MIRRORS**

DESCRIPTION	CODE
Laryngeal Mirrors	75-600-40



**INSUFFLATOR BULB AND TUBE**

DESCRIPTION	CODE
Insufflator Bulb and tube	75-120-11



**TONGUE DEPRESSOR**

DESCRIPTION	CODE
Tongue Depressor With Bulb	75-600-41
Tongue Depressor With LED	75-600-42



**Accessories  
Diagnostic Set  
Capsule and Bulbs**



**75 INDUSTRIES**  
CHOOSE US, TRUST US

Made From Thermoplastic Material.  
Smooth And Shiny Surface  
Ensure Easy Cleaning.  
Bayonet Lock For Otoscope.

**DISPOSABLE SPECULA**

DESCRIPTION	Code
Disposable Specula Dia 2.5mm, Box of 1000	600-43
Disposable Specula Dia 4.0mm, Box of 1000	600-44
Disposable Specula Dia 2.5mm, Box of 250	600-45
Disposable Specula Dia 4.0mm, Box of 250	600-46
Dia 2.5mm s 4. 0mm in Tube	600-47
Clear Tube For Specula	600-48



Made From Thermoplastic Material.  
Smooth And Shiny Surface  
Ensure Easy Cleaning.  
Autoclavable At 130 C.  
Metallic Bayonet Lock For Otoscope.  
Can Be Disinfected.

**REUSABLE SPECULA**

DESCRIPTION	Code
Reusable Specula No.1 Dia 2.5mm	600-49
Reusable Specula No.2 Dia 3.5mm	600-50
Reusable Specula No.3 Dia 4.5mm	600-51
Reusable Specula No.4 Dia 5.5mm	600-52
Reusable Specula No.1 Dia 9.5mm	600-53



Made From Thermoplastic Material.  
Smooth And Shiny Surface  
Ensure Easy Cleaning.  
Autoclavable At 130 C.  
Metallic Bayonet Lock For Otoscope.  
Can Be Disinfected.

**REUSABLE VETERINARY SPECULA**

DESCRIPTION	Code
REUSABLE VETERINARY SPECULA No.1 Dia 4.0mm	600-54
REUSABLE VETERINARY SPECULA No.2 Dia 5.0mm	600-55
REUSABLE VETERINARY SPECULA No.3 Dia 7.0mm	600-56



For Cheap conventional Laryngoscopes blades  
LUX . 16,000

For conventional Laryngoscopes blades  
LUX . 20,000

For Fiber Optic Laryngoscopes blades  
LUX . 20,000 to 40.000

Best Quality LED  
With Almost 50,000 LUX

LED Capsule Fix  
With Almost 50,000 To 100.000 LUX


Item	Code
Bulb W/Lining	75-01
Bulb W/o Lining	75-02

Item	Code
Bulb W/Lining	75-03
Bulb W/o Lining	75-04

Item	Code
LED 20,000L	75-05
LED 40,000L	75-06

Item	Code
LED 35,000L	75-07
LED 50,000L	75-08

Item	Code
LED 50,000L	75-09
LED 100,000L	75-10



For Mini And Metal Body Ophthalmoscope

For 4 Aperture Ophthalmoscope

Xenon

Halogen

Standard Capsule  
With Replaceable Bulb s LED

Item	Code
Bulb	75-11
LED	75-12

Item	Code
Bulb	75-13
LED	75-14

Item	Code
Xenon	75-15

Item	Code
Halogen 2.5 V	75-16
Halogen 3.5 V	75-160

Item	Code
With Bulb	75-17
With LED	75-18



### **31) Diagnostic Hammers**