



- 1) Macintosh Blades - Standard**
- 2) Miller Blades - Standard**
- 3) Macintosh Blades - FIBEROPTIC**
- 4) Macintosh Blades - FIBEROPTIC Fix tube**
- 5) Miller Blades - FIBEROPTIC**
- 6) Miller Blades - FIBEROPTIC Fix tube**
- 7) MEGA Blades - FIBEROPTIC**

## Macintosh Blade Fiber Optic illumination

### Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Available In 3mm, 4mm, 5mm Fiber Optic illumination
- ▶ Inserted Fiber Optic Technology For Advanced Performance
- ▶ Individual Micro-fibers For Improved Light Transmission
- ▶ And Longer Life. Smooth Shape Avoids Sharp Edges & Cavities For Easy Cleaning And Sterilizing.
- ▶ Extra-strong Removable Fiber Bundles To Ensure The
- ▶ Best Possible Light Intensity And Durability.
- ▶ The Clear Design Of The Macintosh F.o. Blades Simplifies Cleaning And Autoclaving. Fiber Optic Pipe Is Removable
- ▶ It Can Be Remove or Replace. Require Fiber Optic Handle.

### Size Measurement

Code	Blades	B/No	Overall length
75-101-01	Premature	00	68 MM
75-101-02	Neonate	0	80 MM
75-101-03	Infant	1	95 MM
75-101-04	Child	2	105 MM
75-101-05	Medium Adult	3	135 MM
75-101-06	Large Adult	4	155 MM
75-101-07	X-Large	5	175 MM

### F.O MAC BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	21-200-30
Set of 3 Blades Size 1,2,3 and Adult Handle	21-200-31
Set of 4 Blades Size 0, 1,2,3 and Adult Handle	21-200-40
Set of 4 Blades Size 1, 2,3,4 and Adult Handle	21-200-41



Half Weld Pipe

Full Weld Pipe

## Macintosh Blade

### Conventional illumination

### Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Available With LED And Standard Bulb Illumination
- ▶ Innovative design for simple incubation and good visibility Autoclave Able
- ▶ Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux
- ▶ These Blades Need No Maintenance — Ever But If Problem With Bulb Then It Can Be Replace.
- ▶ Quick-to-remove Lamp At The Front Of The Blade.
- ▶ The Rounded Edges Of The Blade Reduce The Risk Of Injury. Require Conventional Handle

### Size Measurement

Code	Blades	B/No	Overall length
75-100-00	Premature	00	68 MM
75-100-01	Neonate	0	80 MM
75-100-02	Infant	1	95 MM
75-100-03	Child	2	105 MM
75-100-04	Medium Adult	3	135 MM
75-100-05	Large Adult	4	155 MM
75-100-06	X-Large	5	175 MM

### MACINTOSH CONVENTIONAL SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	75-200-30
Set of 3 Blades Size 1, 2, 3 and Adult Handle	75-200-31
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	75-200-40
Set of 4 Blades Size 1, 2, 3,4 and Adult Handle	75-200-41



## Miller Blade Conventional illumination

Half Weld Pipe

Full Weld Pipe

### Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel.
- ▶ All Blades Compatible With ISO 7376 Standard.
- ▶ Available With Led And Standard Bulb Illumination.
- ▶ Innovative design of the miller blade provides good visibility of the epiglottis & vocal chords. Autoclave Able.
- ▶ Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux.
- ▶ These Blades Need No Maintenance — Ever But If Problem With Bulb Then It Can Be Replace.
- ▶ Compatible With All Conventional Handles.
- ▶ Quick-to-remove Lamp At The Front Of The Blade.
- ▶ The Rounded Edges Of The Blade Reduce The Risk Of Injury. Require Conventional Handle.

### Size Measurement

Code	Blades	B/No	Overall length
75-102-01	Premature	00	65 MM
75-102-02	Neonate	0	78 MM
75-102-03	Infant	1	102 MM
75-102-04	Child	2	154 MM
75-102-05	Medium Adult	3	195 MM
75-102-06	Large Adult	4	205 MM

### For Miller BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 00,0,1 and Handle	22-200-30
Set of 3 Blades Size 0,1,2 and Handle	22-200-31
Set of 4 Blades Size 00,0,1,2 and Handle	22-200-40
Set of 4 Blades Size 0,1, 2, 3 and Handle	22-200-41



## Miller Blade Fiber Optic illumination

### Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Available In 3mm, 4mm, 5mm Fiber Optic illumination
- ▶ Fiber Optic Pipe Is Removable It Can Be Remove Or Replace.
- ▶ Inserted Fiber Optic Technology For Advanced Performance
- ▶ Individual Micro-fibers For Improved Light Transmission
- ▶ And Longer Life. Smooth Shape Avoids Sharp Edges & Cavities For Easy Cleaning And Sterilizing.
- ▶ The Flat Design Of The F.O. Miller Blade Improves Visibility Considerably And Simplifies Incubation.
- ▶ This Blade Reduces The Risk Of Trauma To The Patient And
- ▶ Provides Excellent Visualization. Require Fiber Optic Handle.

### Size Measurement

Code	Blades	B/No	Overall length
75-103-04	Premature	00	65 MM
75-103-05	Neonate	0	78 MM
75-103-06	Infant	1	102 MM
75-103-07	Child	2	154 MM
75-103-08	Medium Adult	3	195 MM
75-103-09	Large Adult	4	205 MM

### F.O Miller BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 00,0,1 and Handle	23-200-30
Set of 3 Blades Size 0,1,2 and Handle	23-200-31
Set of 4 Blades Size 00,0,1,2 and Handle	23-200-40
Set of 4 Blades Size 0,1, 2, 3 and Handle	23-200-41



## MEGA Blade

### Fiber Optic illumination



### Specification:

- ▶ Premium Quality Blades
- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Inserted Fiber Optic Technology For Advanced Performance Up To 6500 Individual Micro-fibers For Improved Light Transmission And Longer Life These Blades Need No Maintenance - Ever. 5 Year Warranty
- ▶ Built With An Integral, Extra-large Fiber Optic Bundle Providing Cooler, White Light.
- ▶ Smooth Shape Avoids Sharp Edges And Cavities For Easy Cleaning And Sterilizing.
- ▶ There Are No Connecting Screws To Get Lost, No External Light Pathways To Get Wrong.
- ▶ Integral Design Enables Superior Cleaning By Elimination Of Cracks And Crevices.
- ▶ Both Ends Of The Blades Are Protected For Life.
- ▶ Require Fiber Optic Handle.

### Size Measurement

Code	Blades	B/No	Overall length
75-104-04	Neonate	0	80 MM
75-104-05	Infant	1	92 MM
75-104-06	Child	2	115 MM
75-104-07	Medium Adult	3	135 MM
75-104-08	Large Adult	4	155 MM
75-104-09	X-Large	5	175 MM

### F.O MEGABLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 00, 0, 1 and Adult Handle	24-200-30
Set of 3 Blades Size 1,2,3 and Adult Handle	24-200-31
Set of 4 Blades Size 0, 1,2,3 and Adult Handle	24-200-40
Set of 4 Blades Size 1, 2,3,4 and Adult Handle	24-200-41

