



75 INDUSTRIES

CHOOSE US, TRUST US

DIAGNOSTIC INSTRUMENTS CATALOGUE

www.azniind.com



ABOUT US



75 INDUSTRIES is known for its quality For Surgical instruments and Laryngoscopes / ENT devices and is an operated business that has been dedicated to providing healthcare professionals with the products they need to perform their art since a long. Our high quality, precision crafted instruments have been an essential product line to our worthy customers. We carry an exclusive range of Surgical Instruments, Surgical Instruments and Healthcare Equipments.

We are committed to providing a product line that guarantees quality and extraordinary value to the end user as well as a level of service, satisfaction, and cooperation that has become a hallmark of **75 INDUSTRIES**. We have gained a reputation for offering unbeatable prices on branded items as well as outstanding customer service on all levels. Every single instrument is manufactured with great care. We are working hard for maintaining higher standard. Advanced techniques are used in manufacturing and finishing departments in order to make best products.

75 INDUSTRIES is offering you huge selection of our products and much more for a fast growing number of customers. We make it easy for you to browse our products, place orders & benefit from our quality services

We assure you quality instruments, most competitive prices and shortest delivery time.

ABOUT LARYNGOSCOPES

What is it ?

A laryngoscope is an instrument used to examine the larynx. The classic laryngoscope is available with a variety of blade attachments and with/without a Fiber Optic light source. The trend regarding blades has/is moving toward the disposable version.

Function:

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube.

Main Feature:

The Miller Laryngoscope Blade . . . The most popular of the straight types. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-intubate patients.

The Macintosh Laryngoscope Blade . . . The Macintosh is the predominate model of the curved blade types. During laryngoscopy, the tip or beak of the blade is compressed into the angle formed by the base of the tongue and the epiglottis, indirectly raising the epiglottis. Some sizes are available flangeless. Also available is a die-stamped or one-piece construction blade with a slightly expanded flange.

CONTENT PAGE



- 1) Macintosh Blades - Standard
- 2) Miller Blades - Standard
- 3) Macintosh Blades - FIBEROPTIC
- 4) Macintosh Blades - FIBEROPTIC Fix tube
- 5) Miller Blades - FIBEROPTIC
- 6) Miller Blades - FIBEROPTIC Fix tube
- 7) MEGA Blades - FIBEROPTIC



- 19) Oto and Ophthalmoscope Re-usable
- 20) Oscope Sets
- 21) Oto and Ophthalmoscope Sets
- 22) Dermatoscopes and ENT set
- 23) Diagnostic Sets
- 24) Mini Oscopes
- 25) Mini Ophthalmoscopes
- 26) Pocket Style Ophthalmoscopes
- 27) Pocket Style Dermatoscopes



- 31) Tuning Forks



- 8) McCoy MEGA Blades - FIBEROPTIC
- 9) McCoy Macintosh Blades - Standard
- 10) McCoy Macintosh Blades - FIBEROPTIC
- 11) Miller Optima CLX Blades - FIBEROPTIC
- 12) Other Blades
- 13) Disposable Blades - Standard
- 14) Disposable Blades - FIBEROPTIC



- 28) Spare accessories Diagnostic set
- 29) Bulbs and Spatulas



- 32) Rigid Laryngoscopes & Microsurgery



- 15) Standard and F.O. Handles
- 16) Standard and F.O. Handles - New Pattern
- 17) Standard and F.O. Handles - Rechargeable
- 18) Charger and Batteries



- 30) Diagnostic Hammers



- 33) Magill Forceps



- 34) Case and Packing



**FIBER OPTIC
GREEN SYSTEM
REUSABLE**

All blades in green system crafted from finest non-magnet stainless steel minimize the risk of contamination.

Cool and brighter light with no chance of tissue burning during long intubation procedures.

Superior quality of glass fiber optic gives maximum bright illumination to view the larynx of human body.

Removable light guide can be autoclave separately, which makes long life of blade.

Autoclave able at 134c°. 5 min, for 2000-3000 cycles at maximum light output.

All blades and handles are made to comply with ISO 7376, allow compatibility with existing blades and handles.

Complete range of handles also equipped with LED lamps.

LED technology has longer additional life time with 50,000hrs and continuous 2hrs operating time.

Brighter than Xenon lamp with 12,000 — 50,000 LUX transmission light at the tip of the blade.

STANDARD CONVENTIONAL REUSABLE BLADE



We have accommodated the whole range of conventional lamp blades into LED choice also, which protects the patient from accidental burns during lengthy procedures.

Standard Vacuum Lamp:

Reliable vacuum standard lamp of 2.5V for better illumination.

Conventional color temperature of 3000k @ 2000 LUX.

2hrs operating time with 20hrs life of lamp.

LED Lamp:

LED lamp gives continuous cool and bright illumination.

2hrs operating time.

Brighter than standard lamp with 12,000 LUX and also 85,000 LUX light beam.

Available with 6500k color temperature.

Longer additional life time with 50,000 hrs.

Whole range of blades complies with ISO 7376 standard.

Easily replace lamp at the front of the blade.

Autoclaved @ 134° C for 5 minutes.

Crafted from finest stainless steel AISI 303 & 304.

Round and smooth edges reduce the risk of injuries.



- 1) Macintosh Blades - Standard**
- 2) Miller Blades - Standard**
- 3) Macintosh Blades - FIBEROPTIC**
- 4) Macintosh Blades - FIBEROPTIC Fix tube**
- 5) Miller Blades - FIBEROPTIC**
- 6) Miller Blades - FIBEROPTIC Fix tube**
- 7) MEGA Blades - FIBEROPTIC**

Macintosh Blade Fiber Optic illumination

Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Available In 3mm, 4mm, 5mm Fiber Optic illumination
- ▶ Inserted Fiber Optic Technology For Advanced Performance
- ▶ Individual Micro-fibers For Improved Light Transmission
- ▶ And Longer Life. Smooth Shape Avoids Sharp Edges & Cavities For Easy Cleaning And Sterilizing.
- ▶ Extra-strong Removable Fiber Bundles To Ensure The Best Possible Light Intensity And Durability.
- ▶ The Clear Design Of The Macintosh F.o. Blades Simplifies Cleaning And Autoclaving. Fiber Optic Pipe Is Removable
- ▶ It Can Be Remove or Replace. Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-101-01	Premature	00	68 MM
75-101-02	Neonate	0	80 MM
75-101-03	Infant	1	95 MM
75-101-04	Child	2	105 MM
75-101-05	Medium Adult	3	135 MM
75-101-06	Large Adult	4	155 MM
75-101-07	X-Large	5	175 MM

F.O MAC BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	21-200-30
Set of 3 Blades Size 1,2,3 and Adult Handle	21-200-31
Set of 4 Blades Size 0, 1,2,3 and Adult Handle	21-200-40
Set of 4 Blades Size 1, 2,3,4 and Adult Handle	21-200-41



Half Weld Pipe

Full Weld Pipe

Macintosh Blade

Conventional illumination

Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Available With LED And Standard Bulb Illumination
- ▶ Innovative design for simple incubation and good visibility Autoclave Able
- ▶ Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux
- ▶ These Blades Need No Maintenance — Ever But If Problem With Bulb Then It Can Be Replace.
- ▶ Quick-to-remove Lamp At The Front Of The Blade.
- ▶ The Rounded Edges Of The Blade Reduce The Risk Of Injury. Require Conventional Handle

Size Measurement

Code	Blades	B/No	Overall length
75-100-00	Premature	00	68 MM
75-100-01	Neonate	0	80 MM
75-100-02	Infant	1	95 MM
75-100-03	Child	2	105 MM
75-100-04	Medium Adult	3	135 MM
75-100-05	Large Adult	4	155 MM
75-100-06	X-Large	5	175 MM

MACINTOSH CONVENTIONAL SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	75-200-30
Set of 3 Blades Size 1, 2, 3 and Adult Handle	75-200-31
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	75-200-40
Set of 4 Blades Size 1, 2, 3,4 and Adult Handle	75-200-41



Miller Blade Conventional illumination

Half Weld Pipe

Full Weld Pipe

Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel.
- ▶ All Blades Compatible With ISO 7376 Standard.
- ▶ Available With Led And Standard Bulb Illumination.
- ▶ Innovative design of the miller blade provides good visibility of the epiglottis & vocal chords. Autoclave Able.
- ▶ Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux.
- ▶ These Blades Need No Maintenance — Ever But If Problem With Bulb Then It Can Be Replace.
- ▶ Compatible With All Conventional Handles.
- ▶ Quick-to-remove Lamp At The Front Of The Blade.
- ▶ The Rounded Edges Of The Blade Reduce The Risk Of Injury. Require Conventional Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-102-01	Premature	00	65 MM
75-102-02	Neonate	0	78 MM
75-102-03	Infant	1	102 MM
75-102-04	Child	2	154 MM
75-102-05	Medium Adult	3	195 MM
75-102-06	Large Adult	4	205 MM

For Miller BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 00,0,1 and Handle	22-200-30
Set of 3 Blades Size 0,1,2 and Handle	22-200-31
Set of 4 Blades Size 00,0,1,2 and Handle	22-200-40
Set of 4 Blades Size 0,1, 2, 3 and Handle	22-200-41



Miller Blade Fiber Optic illumination

Specification:

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Available In 3mm, 4mm, 5mm Fiber Optic illumination
- ▶ Fiber Optic Pipe Is Removable It Can Be Remove Or Replace.
- ▶ Inserted Fiber Optic Technology For Advanced Performance
- ▶ Individual Micro-fibers For Improved Light Transmission
- ▶ And Longer Life. Smooth Shape Avoids Sharp Edges & Cavities For Easy Cleaning And Sterilizing.
- ▶ The Flat Design Of The F.O. Miller Blade Improves Visibility Considerably And Simplifies Incubation.
- ▶ This Blade Reduces The Risk Of Trauma To The Patient And
- ▶ Provides Excellent Visualization. Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-103-04	Premature	00	65 MM
75-103-05	Neonate	0	78 MM
75-103-06	Infant	1	102 MM
75-103-07	Child	2	154 MM
75-103-08	Medium Adult	3	195 MM
75-103-09	Large Adult	4	205 MM

F.O Miller BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 00,0,1 and Handle	23-200-30
Set of 3 Blades Size 0,1,2 and Handle	23-200-31
Set of 4 Blades Size 00,0,1,2 and Handle	23-200-40
Set of 4 Blades Size 0,1, 2, 3 and Handle	23-200-41





MEGA Blade

Fiber Optic illumination

Specification:

- ▶ Premium Quality Blades
- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Inserted Fiber Optic Technology For Advanced Performance Up To 6500 Individual Micro-fibers For Improved Light Transmission And Longer Life These Blades Need No Maintenance - Ever. 5 Year Warranty
- ▶ Built With An Integral, Extra-large Fiber Optic Bundle Providing Cooler, White Light.
- ▶ Smooth Shape Avoids Sharp Edges And Cavities For Easy Cleaning And Sterilizing.
- ▶ There Are No Connecting Screws To Get Lost, No External Light Pathways To Get Wrong.
- ▶ Integral Design Enables Superior Cleaning By Elimination Of Cracks And Crevices.
- ▶ Both Ends Of The Blades Are Protected For Life.
- ▶ Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-104-04	Neonate	0	80 MM
75-104-05	Infant	1	92 MM
75-104-06	Child	2	115 MM
75-104-07	Medium Adult	3	135 MM
75-104-08	Large Adult	4	155 MM
75-104-09	X-Large	5	175 MM

F.O MEGABLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 00, 0, 1 and Adult Handle	24-200-30
Set of 3 Blades Size 1,2,3 and Adult Handle	24-200-31
Set of 4 Blades Size 0, 1,2,3 and Adult Handle	24-200-40
Set of 4 Blades Size 1, 2,3,4 and Adult Handle	24-200-41





- 8) McCoy MEGA Blades - FIBEROPTIC**
- 9) McCoy Macintosh Blades - Standard**
- 10) McCoy Macintosh Blades - FIBEROPTIC**
- 11) Miller Optima CLX Blades - FIBEROPTIC**
- 12) Other Blades**
- 13) Disposable Blades - Standard**
- 14) Disposable Blades - FIBEROPTIC**

FlexTip 70°



McCOY - MEGA Blade

Fiber Optic illumination



Specification:

The McCoy Is A Macintosh Laryngoscope Blade With A Tip Which Is Adjustable Through 70°. A Lever Controls The Tip Angle During Incubation To Lift The Epiglottis, Giving A Clear View Of The Vocal Cords. Contact With The Upper Row Of Teeth And Unnecessary Pressure On The Tongue And Soft Tissue Are Avoided.

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Inserted Fiber Optic Technology For Advanced Performance Up To 6500 Individual Micro-fibers For Improved Light Transmission And Longer Life.
- ▶ These Blades Need No Maintenance — Ever.
- ▶ Built With An Integral, Extra-large Fiber Optic Bundle Providing Cooler, White Light.
- ▶ Smooth Shape Avoids Sharp Edges And Cavities For Easy Cleaning And Sterilizing..
- ▶ Integral Design Enables Superior Cleaning By Elimination Of Cracks And Crevices.
- ▶ Both Ends Of The Blades Are Protected For Life.
- ▶ Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-106-04	Child	2	102 MM
75-106-05	Medium Adult	3	130 MM
75-106-06	Large Adult	4	155 MM

Fo MCCOY BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 2 Blades Size 2 and Adult Handle	75-420-53
Set of 2 Blades Size 3 and Adult Handle	75-420-54
Set of 2 Blades Size2, 3 and Adult Handle	75-420-55
Set of 4 Blades Size2, 3 4 and Child Handle	75-420-56



FlexTip 70°



**McCoy - Blade
Standard illumination**

Specification:

The McCoy Is A Macintosh Laryngoscope Blade With A Tip Which Is Adjustable Through 70°. A Lever Controls The Tip Angle During Incubation To Lift The Epiglottis, Giving A Clear View Of The Vocal Cords. Contact With The Upper Row Of Teeth And Unnecessary Pressure On The Tongue And Soft Tissue Are Avoided.

- ▶ Manufactured From Aisi 304 Stainless Steel.
- ▶ All Blades Compatible With Iso 7376 Standard.
- ▶ Available With Led And Standard Bulb Illumination.
- ▶ Innovative Design For Simple Incubation And Good Visibility.
- ▶ Autoclave Able.
- ▶ Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux.
- ▶ These Blades Need No Maintenance — Ever But If Problem With Bulb Then It Can Be Replace.
- ▶ Quick-to-remove Lamp At The Front Of The Blade.
- ▶ The Rounded Edges Of The Blade Reduce The Risk Of Injury.
- ▶ Require Conventional Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-107-04	Child	2	102 MM
75-107-05	Medium Adult	3	130 MM
75-107-06	Large Adult	4	155 MM

FO MCCOY BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 2 Blades Size 2 and Adult Handle	75-430-53
Set of 2 Blades Size 3 and Adult Handle	75-430-54
Set of 2 Blades Size 2, 3 and Adult Handle	75-430-55
Set of 4 Blades Size 2, 3 4 and Child Handle	75-430-56



FlexTip 70°



**Also
FIX Tube**

McCoy - Blade

Fiber Optic illumination



Specification:

The McCoy Is A Macintosh Laryngoscope Blade With A Tip Which Is Adjustable Through 70°. A Lever Controls The Tip Angle During Incubation To Lift The Epiglottis, Giving A Clear View Of The Vocal Cords. Contact With The Upper Row Of Teeth And Unnecessary Pressure On The Tongue And Soft Tissue Are Avoided.

- ▶ Manufactured From Aisi 304 Stainless Steel
- ▶ All Blades Compatible With Iso 7376 Standard
- ▶ Autoclave Able.
- ▶ Available In 3mm, 4mm, 5mm Fiber Optic Illumination
- ▶ Inserted Fiber Optic Technology For Advanced Performance
- ▶ Individual Micro-fibers For Improved Light Transmission.
- ▶ Smooth Shape Avoids Sharp Edges And Cavities For Easy Cleaning And Sterilizing.
- ▶ Extra-strong Fiber Bundles To Ensure The Best Possible Light Intensity And Durability.
- ▶ The Clear Design Of The Macintosh F.O. Blades Simplifies Cleaning And Autoclaving
- ▶ Fix Tube - Fiber Optic Pipe Is Non-removable It Cant Be Remove.
- ▶ Also Available In Standard Fiber Optic Illumination (removable pipe). Require Fiber Optic Handle.

Size Measurement

Code	Blades	B/No	Overall length
75-105-04	Child	2	105 MM
75-105-05	Medium Adult	3	130 MM
75-105-06	Large Adult	4	155 MM

Fo MCCOY BLADE SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 2 Blades Size 2 and Adult Handle	75-440-53
Set of 2 Blades Size 3 and Adult Handle	75-440-54
Set of 2 Blades Size2, 3 and Adult Handle	75-440-55
Set of 4 Blades Size2, 3 4 and Child Handle	75-440-56



Disposable Blades Standard illumination

Specification:

- ▶ Single-use, sterile disposable blades help eliminate the time, reprocessing steps and costs related to disinfecting reusable blades.
- ▶ Packaged individually in a clear pouch making it easy to identify blade size and type.
- ▶ Reduced risk of cross-contamination
- ▶ Stainless steel body
- ▶ Macintosh sizes 0-4
- ▶ Miller sizes 00-4
- ▶ Interchangeable with all standard handles



Miller Size Measurement

Code	Blades	B/No	Overall length
75-200-93	Premature	00	68 MM
75-200-94	Neonate	0	80 MM
75-200-95	Infant	1	95 MM
75-200-96	Child	2	105 MM
75-200-97	Medium Adult	3	135 MM
75-200-98	Large Adult	4	155 MM

Mac Size Measurement

Code	Blades	B/No	Overall length
75-200-99	Neonate	0	80 MM
75-200-00	Infant	1	95 MM
75-200-01	Child	2	105 MM
75-200-02	Medium Adult	3	135 MM
75-200-03	Large Adult	4	155 MM

MACINTOSH CONVENTIONAL SETS

Standard Sets (Custom Sets Available)

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	75-450-57
Set of 3 Blades Size 1, 2, 3 and Adult Handle	75-450-58
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	75-450-59
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	75-450-60

Disposable Blades Standard illumination

Specification:

- ▶ Single-use, sterile disposable blades help eliminate the time, reprocessing steps and costs related to disinfecting reusable blades.
- ▶ Packaged individually in a clear pouch making it easy to identify blade size and type.
- ▶ Reduced risk of cross-contamination
- ▶ Stainless steel body
- ▶ Macintosh sizes 0-4
- ▶ Miller sizes 00-4
- ▶ Interchangeable with all green-band fiber optic handles



Miller Size Measurement

Code	Blades	B/No	Overall length
75-800-05	Premature	00	68 MM
75-800-06	Neonate	0	80 MM
75-800-07	Infant	1	95 MM
75-800-08	Child	2	105 MM
75-800-09	Medium Adult	3	135 MM
75-800-10	Large Adult	4	155 MM


Mac Size Measurement

Code	Blades	B/No	Overall length
75-800-11	Neonate	0	80 MM
75-800-12	Infant	1	95 MM
75-800-13	Child	2	105 MM
75-800-14	Medium Adult	3	135 MM
75-800-15	Large Adult	4	155 MM

MACINTOSH CONVENTIONAL SETS

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	75-460-57
Set of 3 Blades Size 1, 2, 3 and Adult Handle	75-460-58
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	75-460-59
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	75-460-60



- 15) Standard and F.O. Handles**
- 16) Standard and F.O. Handles - New Pattern**
- 17) Standard and F.O. Handles Rechargeable**
- 18) Charger and Batteries**

Replacements Bulbs

Bulb	Battery	Code
Xenon	2.5V	75-15
LED	3.0V	75-07
Halogen	2.5V	75-16
Halogen	3.5V	75-160



**Handles
Fiber Optic Illumination**

Specification:

- ▶ Chrome Plated Brass Fiber Optic Handle compatible with all "Green System" fiber optic blades.
- ▶ Require 2 AA size Batteries or 3.5 V Rechargeable Battery.
- ▶ Crafted From High quality Chrome plated Brass.
- ▶ Autoclaveable

Battery Handle Fiberoptic Xenon, Brass Chrome Plated

Handle	Overall length	Battery	Code
Infant	122 MM	1 x AA Size	75-300-17
Stubby	125 MM	1 x C Size	75-300-18
Pediatric	160 MM	2 x AA Size	75-300-19
Medium	160 MM	2x C size	75-300-20
Large	175 MM	2x D size	75-300-21

Battery Handle Fiberoptic LED, Brass Chrome Plated

Handle	Overall length	Battery	Code
Infant	122 MM	1 x AA Size	75-300-22
Stubby	125 MM	1 x C Size	75-300-23
Pediatric	160 MM	2 x AA size	75-300-24
Medium	150 MM	2x C size	75-300-25
Large	175 MM	2x D size	75-300-26

**Handles
Conventional Illumination**

Battery Handle Conventional, Brass Chrome Plated


Handle	Overall length	Battery	Code
Infant	110 MM	1 x AA Size	75-300-27
Stubby	115 MM	1 x C Size	75-300-28
Pediatric	145 MM	2 x AA size	75-300-29
Medium	150 MM	2 x C size	75-300-30
Large	165 MM	2 x D size	75-300-31



Handles
Fiber Optic Illumination

Handle	Overall length	Code
Pediatric	160 MM	75-300-32
Medium	150 MM	75-300-33
Medium With DC	150 MM	75-300-33-1

Rechargeable Handle
Bottom For Desktop Charger



DC Rechargeable Handle
Bottom



Handles
Conventional Illumination

Handle	Overall length	Code
Pediatric	145 MM	75-300-34
Medium	150 MM	75-300-35
Medium With DC	150 MM	75-300-35-1

Handles
Oto & Ophthalmoscope

Handle	Overall length	Code
Aluminum Handle	145 MM	75-300-36
Medium	150 MM	75-300-37

- ▶ Chrome Plated Brass Fiber Optic Handle compatible with all "Green System" fiber optic blades.
- ▶ 3.5 V Rechargeable Battery.
- ▶ Crafted From High quality Chrome plated Brass.
- ▶ Autoclaveable



Desk Charger Twin Port For Mega and Baby Handle

Specification:

- ▶ Manufactured From Plastic Material
- ▶ Big Handles requires 3.5 V size Rechargeable batteries
- ▶ Baby Handles requires 3 V size Rechargeable batteries

Purpose And Features

- ▶ The DUKE is only intended to charge rechargeable battery handles.
- ▶ Power-on indicator. Charging state indicator.
- ▶ When take the handle away the charger indicator
- ▶ Automatically turns off
- ▶ 3. 7V Lithium battery.
- ▶ 2X charging ports.
- ▶ Circuit short, reverse-voltage, over-voltage protection.
- ▶ Low battery indicator.

Use This For Child Handle
Converter for Rechargeable Baby Handle
Charge with Desk Charger



300-38-DUKE
With 9-12 V DC charger
Every Port Provide 4.5 V OUTPUT

Li-ion Batteries



LI-ION BATTERY

For Big Handle 4.0 V
Can Only Charge With Desk
Charger



DC LI-ION BATTERY

For Diagnostics Set & Big Handle
Can Only Charge With Normal
Wall Charger Of 5-6 V
Also Available For Standard And
Fiber Optic Handle



LI-ION BATTERY

For Baby Handle 3.0 V For
Standard And Fiber Optic
Can Only Charge With Desk
Charger

DESCRIPTION	Code
Adult & Otoscope Handle 2.5V	75-300-39
Adult & Otoscope Handle 3.5V	75-300-40

DESCRIPTION	Code
Adult & Otoscope Handle 2.5V	75-300-41
Adult & Otoscope Handle 3.5V	75-300-42

DESCRIPTION	Code
Adult & Small Handle 2.5V	75-300-43
Adult & Small Handle 3.5V	75-300-44



- 19) Oto and Ophthalmoscope Re-usable**
- 20) Otoscope Sets**
- 21) Oto and Ophthalmoscope Sets**
- 22) Dermatoscopes and ENT set**
- 23) Diagnostic Sets**
- 24) Mini Otoscopes**
- 25) Mini Ophthalmoscopes**
- 26) Pocket Style Ophthalmoscopes**
- 27) Pocket Style Dermatoscopes**

Otosopes and Ophthalmoscopes

Good Quality and Reusable



75-200-04
Brass Otoscope Fiber Optic
Brass Otoscope Conventional



75-200-05
ENT Set Complete
Pouch Packing



75-200-06
Brass Ophthalmoscope

Specification:

- ▶ Manufactured From Brass Chrome plated
- ▶ 3x Magnifying glass
- ▶ Available in Halogen and LED Illumination
- ▶ Bulbs 2.5V
- ▶ Handles requires 2 x C size batteries
- ▶ Hard Case and Soft Pouch available
- ▶ Fiber Optic Light Source

Fittings



Pin Fitting



Bayonet Fitting

Otoscopes And Ophthalmoscopes STANDARD SETS

BASIC OTOSCOPE SET Bayonet Fitting

Complete With Otoscope, Battery Handle, 1 Set Of 5 Reusable Specula Dia 2.5, 3.5, 4.5, 5.5, 9.0mm In a Pouch

CODE	DESCRIPTION
75-500-72	Completes Set With Spare Lamp
75-500-73	Complete Set Without Spare Lamp
75-500-74	Complete Set With Spare Lamp And Rechargeable Battery
75-500-75	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-76	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger

BASIC OPHTHALMOSCOPE SET

Complete With OPHTHALMOSCOPE, Battery Handle, In a Pouch

CODE	DESCRIPTION
75-500-77	Completes Set With Spare Lamp
75-500-78	Complete Set Without Spare Lamp
75-500-79	Complete Set With Spare Lamp And Rechargeable Battery
75-500-80	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-81	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger

STANDARD OTOSCOPE SET Pin Fitting

Complete With Otoscope, Battery Handle, 1 Set Of 5 Reusable Specula Dia 2.5, 3.5, 4.5, 5.5, 9.0mm In a Pouch

CODE	DESCRIPTION
75-500-82	Completes Set With Spare Lamp
75-500-83	Complete Set Without Spare Lamp
75-500-84	Complete Set With Spare Lamp And Rechargeable Battery
75-500-85	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-86	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger

BASIC NASAL OTOSCOPE SET Pin Fitting

Complete With Otoscope, Battery Handle, 1 Nasal Speculum, In A Pouch
Also Available In Plastic Box

CODE	DESCRIPTION
75-500-87	Completes Set With Spare Lamp
75-500-88	Complete Set Without Spare Lamp
75-500-89	Complete Set With Spare Lamp And Rechargeable Battery
75-500-90	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-91	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



Otosopes And Ophthalmoscopes STANDARD SETS

BASIC OPHTHALMOSCOPE SET

Complete With Ophthalmoscope, Battery Handle, In A Plastic Box

CODE	DESCRIPTION
75-500-92	Completes Set With Spare Lamp
75-200-06	Complete Set Without Spare Lamp
75-500-94	Complete Set With Spare Lamp And Rechargeable Battery
75-500-95	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-96	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



STANDARD OTOSCOPE AND OPHTHALMOSCOPE SET

Complete With Otoscope ,ophthalmoscope, Battery Handle, In A Plastic Box

CODE	DESCRIPTION
75-500-97	Completes Set With Spare Lamp
75-500-98	Complete Set Without Spare Lamp
75-500-99	Complete Set With Spare Lamp And Rechargeable Battery
75-500-00	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-01	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



STANDARD OTOSCOPE SET

Complete With Otoscope, Battery Handle, 1 Set Of 3 Reusable Specula Dia 4.5, 5.5, 9.0mm In Plastic Box

CODE	DESCRIPTION
75-500-02	Completes Set With Spare Lamp
75-200-04	Complete Set Without Spare Lamp
75-500-04	Complete Set With Spare Lamp And Rechargeable Battery
75-500-05	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-06	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



pin Fitting

Bayonet Fitting

BASIC DERMATOSCOPES SET

Complete With Dermatoscopes, Spare Mirror, In A Plastic Box

CODE	DESCRIPTION
75-500-07	Completes Set With Spare Lamp
75-500-08	Complete Set Without Spare Lamp
75-500-09	Complete Set With Spare Lamp And Rechargeable Battery
75-500-10	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-11	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



BASIC DERMATOSCOPES SET

Complete With Dermatoscopes, Spare Mirror, In A Pouch

CODE	DESCRIPTION
75-500-12	Completes Set With Spare Lamp
75-500-13	Complete Set Without Spare Lamp
75-500-14	Complete Set With Spare Lamp And Rechargeable Battery
75-500-15	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-16	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



STANDARD DIAGNOSTIC SET

Complete With Standard Otoscope, Standard Ophthalmoscope, Handle, Nasal speculum, Laryngeal mirrors, Reusable specula, Hard case

CODE	DESCRIPTION
75-500-17	Completes Set With Spare Lamp
75-500-18	Complete Set Without Spare Lamp
75-500-19	Complete Set With Spare Lamp And Rechargeable Battery
75-500-20	Complete Set With Spare Lamp And Rechargeable Battery And Desk Charger
75-500-21	Complete Set With Spare Lamp And Rechargeable Battery And Wall Charger



Diagnostic Set Full Specification

UNIVERSAL DIAGNOSTIC SET

75-2000-04



- Handle (dry Cell Or Rechargeable) X 1
- Otoscope
- Ophthalmoscope
- Laryngeal Mirrors X2
- Reusable Specula X 5
- Tongue Blade Holder
- Metal Tongue Blade Depressor
- Nasal Speculum
- Ear Tip Connector
- Ear Tips
- Replacement Bulb X 2
- Hard Case X 1
- Also Available In Leather Pouch Just Add -1 After Code

RIESTER DIAGNOSTIC SET

75-2000-02



- Standard Otoscope X 1
- Standard Ophthalmoscope X 1
- Handle (dry Cell Or Rechargeable) X 1
- Nasal Speculum
- Laryngeal Mirrors 2x
- Reusable Specula X 5
- Replacement Bulb X 2
- Hard Case X 1
- Also Available In Leather Pouch Just Add -1 After Code

More Details

- Manufactured From High Quality Material
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available

About Sets

Diagnostic Sets Are The Most Commonly Used Tool Sets During A Clinicians Repertoire To Perform A Comprehensive Physical Test. Diagnostic Sets Facilitate Physicians To Perform Quick Exams And Accurately.



75-300-04
 Mini Otoscope Fiber Optic
 Complete Set of Disposable Cannulas



75-300-05
 Mini Otoscope Fiber Optic
 Complete Set of Disposable Cannulas

Specification:

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source



A B C D

Color Options





75- 500-22



75- 500-20

Mini Ophthalmoscope

Parker mini Ophthalmoscope with 2.5v LED bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses - 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification fo adopter setting.

Power Supply: 2 x AA Batteries (not included in set).
Packing: Hard plastic case packing.



75-500-23



75-500-24



75- 500-25



Multi-Action Tongue Depressor

All-Purpose Pocket-Sized Illuminator that is Easy to Use and Maintain.
Multiple Applications.

Tongue-Depressor Holder can be Used with Disposable Depressors (Included) or with Standard Wooden Depressors. Bright Halogen Bulb with More Intensity Compared with Conventional Bulbs. Refined Plastic Compact Body for Long Lasting Use. Attachment Clip with On/Off Button.

Set Includes: (1) Handle, (1) Head with Integrated Tongue-Depressor Holder (3) Disposable Depressors LED Illuminating Probe. Blue Leather Case. Two "AA" Batteries Required.
Batteries not Included.

The ophthalmoscope is a tool used in medical to examine the interior of the eye. Studying the inner parts of the eye, such as the optic disk, retinal blood vessels, retina, choroid and macula, can help diagnose diseases of the eye.



Pocket Ophthalmoscope

75-450

The Pocket size Ophthalmoscope is a lightweight, compact diagnostic instrument that delivers high performance in a slightly smaller package. Quality optics in mini-format is ideal for rounds and house calls. Attachment clip with integrated on/off switch. High-quality plastic handle with shockproof, sturdy, non-slip features. Enhanced LED Technology for perfect illumination and a brilliant image. 4 different apertures. Comprehensive lens range from + 20 D to -20 D.

Power Supply: 2 x AA Batteries (not included in set).
Packing: Hard plastic case packing.



Pocket Ophthalmoscope

75-500

2.5v LED lamp light for true tissue color and consistent long lasting illumination 16 color coded lenses -25 to +25 for excellent resolution (Red for - Black for +) 4 different apertures (small, large, blue, half moon focus)
Rubber eye brow rest to prevent scratching of eye-glasses
illumination lens dial for clear identification of adopter setting.

Power Supply: 2 x AA Batteries (not included in set).
Packing: Hard plastic case packing.



Dermatoscope

75-550

Modern pocket Dermatoscope in a compact design with high-quality optics. 10X magnification large, high-resolution image. Choice of two skin-friendly and sterilizable contact plates: with scale from 0 - 10mm for exact measurement of pigmented skin lesions, without scale. Simple lamp change on the underside of the instruments head. Enhanced LED Technology for more light. Attachment clip with integrated On/Off switch. High-quality 2 color handle.

Power Supply: 2 x AA Batteries (not included in set).
Packing: Hard plastic case packing.



Dermatoscope

75-610

Bayonet locking system for easy engagement and disengagement from the handle.

2.5v LED illumination of the examination field. High resolution achromatic lenses system with 10x magnification. Adjustable focusing Eyepiece for a large correction range. Lens without scale, Soft eyepiece protection.

Available in Polish and Sand finish.

Power Supply: 2 x C Batteries (not included in set).

Packing: Hard plastic case packing.



- 29) Spare accessories Diagnostic set**
- 30) Bulbs and Spatulas**

TONGUE DEPRESSOR

DESCRIPTION	CODE
Tongue Depressor	75-600-21



CURVED ILLUMINATOR

DESCRIPTION	CODE
Bayonet Fitting	75-600-32
Pin Fitting	75-600-33
Re-usable	75-600-34



METAL TONGUE DEPRESSOR

DESCRIPTION	CODE
Tongue Depressor	75-600-22



OPHTHALMOSCOPE HEAD

DESCRIPTION	CODE
Bayonet Fitting	75-600-35
Pin Fitting	75-600-36
Re-usable	75-600-37



OTOSCOPE HEAD

DESCRIPTION	CODE
W/O Ear Spatula	75-600-23
With Ear Spatula	75-600-24
Bayonet Fitting	75-600-25
Screw Fitting	75-600-26
With LED	75-600-27
With Bulb	75-600-28



NASAL SPECULUM

DESCRIPTION	CODE
Nasal Speculum	75-600-38



ADAPTER FOR DISPOSABLE T.D

DESCRIPTION	CODE
Adapter For Disposable T.D	75-600-29



SPREADABLE NASAL SPECULUM

DESCRIPTION	CODE
SPREADABLE NASAL SPECULUM	75-600-39



METAL ADAPTER FOR SPECULA

DESCRIPTION	CODE
Metal Adapter For BY Fitting	75-600-30
Metal Adapter For PN Fitting	75-600-31



LARYNGEAL MIRRORS

DESCRIPTION	CODE
Laryngeal Mirrors	75-600-40



INSUFFLATOR BULB AND TUBE

DESCRIPTION	CODE
Insufflator Bulb and tube	75-120-11



TONGUE DEPRESSOR

DESCRIPTION	CODE
Tongue Depressor With Bulb	75-600-41
Tongue Depressor With LED	75-600-42



**Accessories
Diagnostic Set
Capsule and Bulbs**



75 INDUSTRIES
CHOOSE US, TRUST US

Made From Thermoplastic Material.
Smooth And Shiny Surface
Ensure Easy Cleaning.
Bayonet Lock For Otoscope.

DISPOSABLE SPECULA

DESCRIPTION	Code
Disposable Specula Dia 2.5mm, Box of 1000	600-43
Disposable Specula Dia 4.0mm, Box of 1000	600-44
Disposable Specula Dia 2.5mm, Box of 250	600-45
Disposable Specula Dia 4.0mm, Box of 250	600-46
Dia 2.5mm s 4. 0mm in Tube	600-47
Clear Tube For Specula	600-48



Made From Thermoplastic Material.
Smooth And Shiny Surface
Ensure Easy Cleaning.
Autoclavable At 130 C.
Metallic Bayonet Lock For Otoscope.
Can Be Disinfected.

REUSABLE SPECULA

DESCRIPTION	Code
Reusable Specula No.1 Dia 2.5mm	600-49
Reusable Specula No.2 Dia 3.5mm	600-50
Reusable Specula No.3 Dia 4.5mm	600-51
Reusable Specula No.4 Dia 5.5mm	600-52
Reusable Specula No.1 Dia 9.5mm	600-53



Made From Thermoplastic Material.
Smooth And Shiny Surface
Ensure Easy Cleaning.
Autoclavable At 130 C.
Metallic Bayonet Lock For Otoscope.
Can Be Disinfected.

REUSABLE VETERINARY SPECULA

DESCRIPTION	Code
REUSABLE VETERINARY SPECULA No.1 Dia 4.0mm	600-54
REUSABLE VETERINARY SPECULA No.2 Dia 5.0mm	600-55
REUSABLE VETERINARY SPECULA No.3 Dia 7.0mm	600-56



For Cheap conventional Laryngoscopes blades
LUX . 16,000



For conventional Laryngoscopes blades
LUX . 20,000



For Fiber Optic Laryngoscopes blades
LUX . 20,000 to 40.000



Best Quality LED With Almost 50,000 LUX



LED Capsule Fix With Almost 50,000 To 100.000 LUX


Item	Code
Bulb W/Lining	75-01
Bulb W/o Lining	75-02

Item	Code
Bulb W/Lining	75-03
Bulb W/o Lining	75-04

Item	Code
LED 20,000L	75-05
LED 40,000L	75-06

Item	Code
LED 35,000L	75-07
LED 50,000L	75-08

Item	Code
LED 50,000L	75-09
LED 100,000L	75-10



For Mini And Metal Body Ophthalmoscope

Item	Code
Bulb	75-11
LED	75-12



For 4 Aperture Ophthalmoscope

Item	Code
Bulb	75-13
LED	75-14



Xenon

Item	Code
Xenon	75-15



Halogen

Item	Code
Halogen 2.5 V	75-16
Halogen 3.5 V	75-160



Standard Capsule With Replaceable Bulb s LED

Item	Code
With Bulb	75-17
With LED	75-18



31) Diagnostic Hammers

TAYLOR HAMMER

DESCRIPTION	Code
Taylor Hammer S.S, 17cm	75-700-57
Taylor Hammer S.S, 20cm	75-700-58

The most popular Hammer
Stainless steel Handle
Latex Free
Available in Two Size



TROMNER HAMMER

DESCRIPTION	Code
Tromner Hammer, 175 Gm, S.S, 24cm	75-700-59


S.S. Handle With 2 Sided Synthetic
Rubber Head
LATEX Free
Available In Size 24cm



BERLINER HAMMER

DESCRIPTION	Code
Berliner Hammer, S.S, 17.0cm	75-700-60

Stainless steel Handle also with
plastic grip
Latex Free
Available Size 17 cm



DEJERINE HAMMER

DESCRIPTION	Code
Dejerine Hammer, S.S, 20.0cm	75-700-61
Dejerine Hammer with Needle, S.S, 21.0cm	75-700-62

S.S. Handle With 2 Sided Synthetic
Rubber Head
Latex Free
Handle Conceals Screw-in-bush
Needle Contained Within The Head



BUCK HAMMER

DESCRIPTION	Code
Buck Hammer, S. S, 15.0cm	75-700-63
Buck Hammer, S. S, 18.0cm	75-700-64

S.S. Handle With 2 Sided Synthetic
Rubber Head
Latex Free
Handle Conceals Screw-in-bush
Needle Contained Within The Head



BABINSKI HAMMER

DESCRIPTION	Code
BABINSKI HAMMER With Needle, S.S, 22.0cm	75-700-65

Octagonal Stainless steel Handle with
Concealed Needle
Rubber Head Can Only Used in One
Direction
Latex Free



BABINSKI HAMMER

DESCRIPTION	Code
BABINSKI HAMMER With Needle, S.S, 22.0cm	75-700-66

Octagonal Stainless steel Handle with
Concealed Needle
Rubber Head Can Used in Both
Directions Vertical s Horizontal
Latex Free





PROFESSIONAL DIAGNOSTIC
INSTRUMENTS



32) Tuning Forks

Tuning Forks Full Specification



HARTMAN

Hz	Code
128	75-19
256	75-20
512	75-21
1024	75-22
2048	75-23
4096	75-24



WELDON ALLOY

Hz	Code
128	75-25
256	75-26
512	75-27
1024	75-28
2048	75-29
4096	75-30



WELDON ALLOY with Weight

Hz	Code
128	75-31
256	75-32



KONIG

Hz	Code
128	75-33



LUCAE

Hz	Code
256	75-34
512	75-35



RYDEL SEIFFERT

Hz	Code
64	75-36



TAYLOR

Hz	Code
128	75-37

SPECIFICATION:

Manufactured From Aluminum Also Available in Steel.
 Pouch Packing Available For Single Piece.
 Set Available In Hard Case Or Plastic Box Or Wooden Box.
 Custom Set Available.
 All Tuning Forks Available In 6 Different Color.

ABOUT TUNING FORKS

The Fundamental Mode Of Vibration Is The Mode Most Commonly Associated With Tuning Forks, The Two Tines Of The Fork Alternately Move Toward And Away From Each Other, Each Bending Like A Cantilever Beam, Fixed At The Stem And Free At The Other End. This Is A Symmetric Mode, Since The Two Tines Are Mirror Images Of Each Other






33) Rigid Laryngoscopes & Microsurgery

Jackson Laryngoscope
with handle and light carrier, satin finish inside

Size 1 (12.0 cm)
75-8175-1

Size 2 (17.1 cm)
75-8175-2

Size 3 (18.0 cm)
75-8175-3




Jackson Laryngoscope
with handle and light carrier, satin finish inside

Size 1 (12.0 cm)
75-8175-4

Size 2 (17.1 cm)
75-8175-5

Size 3 (18.0 cm)
75-8175-6



KLEINSASSER laryngoscope
anterior commissure, 22 cm
especially for long neck

Size 1 (13.0 cm)
75-8175-7

Size 2 (15.0 cm)
75-8175-8

Size 3 (17.0 cm)
75-8175-9

Size 4 (22.0 cm)
75-8175-10



Esophogoscope With
Fiber Optic Light Carrier


25 cm
35 cm
45 cm
52 cm

Esophogoscope

Bronchoscope

Negus Aspirating
Bronchoscope

Bronchoscope
Adolescent
Small 25 cm
Medium 30 cm
Large 40 cm



Rigid laryngoscope



34) Magill Forceps

Magill Forceps Full Specification



AK 900/1
17 cm

AK 900/2
20 cm

AK 900/3
25 cm

About Forceps

A Bent Blunt Instrument Used To Facilitate Nasotracheal Intubation. Angled Forceps Used To Guide A Tracheal Tube Into The Larynx or a Nasogastric Tube Into The Esophagus Under Direct Vision. It Is Also Used To Place Pharyngeal Packs And Remove Foreign Bodies.





35) Case and Packing

Packing Options



Made From Thermoplastic.
For All Kinds Of Mac, Miller And
Diagnostic Sets.
Available In 6 Colors.

Description	Code
Pastic Box	75-110-67



Made From Thermoplastic.
For All Kinds Of Mac, Miller And
Diagnostic Sets.
Available In 6 Colors.

Description	Code
pastic Box	75-110-68



Made From Thermoplastic.
For All Kinds Of Mac, Miller And
Diagnostic Sets.
Only In Black Colors.

Description	Code
pastic Box	75-110-69



Made From Leather.
For All Kinds Of Mac, Miller And
Diagnostic Sets.
Inside Colors Of Your Choice.

Description	Code
Leather Pouch	75-110-70




Specially For Mini Oscope And
Ophthalmoscope. Also Available
For All Kinds Of Handle.
Available In All Color Colors.

Description	Code
Pouch	75-110-71




Made From Jeans.
For All Kinds Of Mac, Miller And
Diagnostic Sets.
Available In All Color Colors.

Description	Code
Bag For Set	75-110-72



Blade Packing Box.

Description	Code
Blade No 00	75-110-73
Blade No 0	75-110-74
Blade No 1	75-110-75



Blade Packing Box.

Description	Code
Blade No 2	75-110-76
Blade No 3	75-110-77
Blade No 4	75-110-78
Blade No 5	75-110-79



Box.Packing
For All kinds of MAC, Miller and
Diagnostic Sets In 110-67 Plastic box

Description	Code
Pastic Box	75-110-80

NOTE
Box Design Change Every Year

NOTE
Add -1 After code For Miller Blades
Available For MAC And Miller Blades Only



75 INDUSTRIES

CHOOSE US, TRUST US

DIAGNOSTIC INSTRUMENTS CATALOGUE

